


## PISTON REPLACEMENT KIT

Both pistons PN 133 213 and PN 124 688 have been obsoleted with a new piston (PN 126 853). When either old-style piston is replaced for the first time, the old-style spider (PN 124 691), shovel rod (PN 124 674), and upper check rod (PN 124 673) must also be replaced. Piston replacement kit PN 126 831 contains all the necessary parts to make this conversion. The new items (1-3) are depicted below. After the old-style parts have been replaced, items 1 through 3 can be ordered separately, or replaced as a group, by reordering the kit.



Ref.	Note	Part No.	Description	Qty.
-		126 831	Kit, Replacement, Piston	1
1		126 853	• Piston	1
2		126 857	• Rod, Upper Check	1
3		126 856	• Spacer, 1.81 O.D. x .93 x .250	1