





25B 16:1 Aluminum Pump

Service Kits and Parts

Refer to 25B Pump Manual 1044132 for complete service and parts information.

For parts and technical support call (800) 433-9319.

Air Motor Kits A - 106274 Seal Kit  Solvent Chamber Fluids 248831 Type S (solvents) 900255 Type K (solvents) 140029 Type T (waterbornes) 244854 Type Q (waterbornes)	Air Valve Kits B - 106103 Seal Kit C - 106314 Crankshaft and Arm Kit 	Hydraulic Section Kits D - 503132 Pressure Ball and Seat Kit E - 244826 Siphon Ball and Seat Kit F - 106107 Pressure and Siphon Kit G - Upper Packing Glands H - Lower Packing Glands J - Hydraulic Service Kits  Type G Gland Kit (Housings not included in kit) 				
Packing Types A - Buna-N / Graphite D - Buna-N / PTFE F - PTFE / Leather G - Polyurethane U - UHMWPE	Glands and Service Kits Upper Packing Gland (G) Lower Packing Gland (H) Hydraulic Service Kit (J)	Type A — — 106106	Type D — — 106112	Type F — — 106292	Type G 503216 503212 106391	Type U — — 335380