

25B 16:1 Stainless Steel Pump

Service Kits and Parts

Refer to 25B Pump Manual 1044133 for complete service and parts information.
For parts and technical support call (800) 433-9319.

Air Motor Kits A - 106274 Seal Kit Solvent Chamber Fluid 248831 Type S (solvents) 900255 Type K (solvents) 140029 Type T (waterbornes) 244854 Type Q (waterbornes) Packing Types A - Buna-N, Graphite D - Buna-N, PTFE F - PTFE/Leather G - Polyurethane U - UHMWPE	Air Valve Kits B - 106103 Seal Kit C - 106314 Crankshaft and Arm Kit Hydraulic Section Kits D - 503385 Pressure Ball and Seat Kit (247238 General Finishing Pump) E - 1047077 Pressure Ball and Seat Kit (149570 CleanSpray® Pump) F - 244823 Siphon Kit G - Upper Packing Glands H - Lower Packing Glands J - Packing Service Kits Type G Gland Kit (Housings not included in kit) Included in Kit Upper 986707 Retaining Ring (not included) 986802 Retaining Ring (not included) Included in Kit Lower	Type G Gland Kit (Housings not included in kit) Included in Kit Upper 986707 Retaining Ring (not included) 986802 Retaining Ring (not included) Included in Kit Lower
	Packing Glands and Kits Type AType DType FType GType U Upper Packing Gland (G)-----248050----- Lower Packing Gland (H)-----248051----- Packing Gland Kits (J)160241416024151602416106392335382	