












A7A Single and Dual Spray Guns

Service Kits and Parts

Refer to A7A Manual 1072640 for complete service and parts information.
For parts and technical support call (800) 433-9319.

Piston and Seals Kit			Ball and Seat/Needle and Seat Kits																																															
A - 106230 Use this kit on all standard pressure (1700 psi) A7A guns. 			B, C 		D 	E 	F 	G, H 		I 	J 	K 																																						
																																																		
			Kits: <table><tr><td>B</td><td>106231</td><td>Ball and Seat</td><td>153150</td></tr><tr><td>C</td><td>320031</td><td>Ball and Seat</td><td>342276, 246789, 337303, 1619165</td></tr><tr><td>D</td><td>106280</td><td>Needle and Seat</td><td>248445, 153160</td></tr><tr><td>E</td><td>242575</td><td>Ball and Seat</td><td>510010</td></tr><tr><td>F</td><td>106234</td><td>Needle and Seat</td><td>245375</td></tr><tr><td>G</td><td>159447</td><td>Ball and Seat</td><td>160062</td></tr><tr><td>H</td><td>337305</td><td>Ball and Seat</td><td>342300, 337302, 337304</td></tr><tr><td>I</td><td>106416</td><td>Needle and Seat</td><td>248937, 149097, 160063</td></tr><tr><td>J</td><td>134972</td><td>Ball and Seat w/Retainer</td><td>160062</td></tr><tr><td>K</td><td>160354</td><td>Needle and Seat w/Retainer</td><td>160937, 149097, 160063</td></tr></table>								B	106231	Ball and Seat	153150	C	320031	Ball and Seat	342276, 246789, 337303, 1619165	D	106280	Needle and Seat	248445, 153160	E	242575	Ball and Seat	510010	F	106234	Needle and Seat	245375	G	159447	Ball and Seat	160062	H	337305	Ball and Seat	342300, 337302, 337304	I	106416	Needle and Seat	248937, 149097, 160063	J	134972	Ball and Seat w/Retainer	160062	K	160354	Needle and Seat w/Retainer	160937, 149097, 160063
			B	106231	Ball and Seat	153150																																												
			C	320031	Ball and Seat	342276, 246789, 337303, 1619165																																												
D	106280	Needle and Seat	248445, 153160																																															
E	242575	Ball and Seat	510010																																															
F	106234	Needle and Seat	245375																																															
G	159447	Ball and Seat	160062																																															
H	337305	Ball and Seat	342300, 337302, 337304																																															
I	106416	Needle and Seat	248937, 149097, 160063																																															
J	134972	Ball and Seat w/Retainer	160062																																															
K	160354	Needle and Seat w/Retainer	160937, 149097, 160063																																															
Used On Guns:																																																		
Packing Cartridges																																																		
	153152	Standard	153150, 510010, 248937 709831, 160062																																															
	153090	Stainless Steel	248445, 342276, 153160 337303, 245375, 342300 149097, 337302, 160063 337304, 1619165																																															
	245778	Corrosion Resistant	246789																																															
	244918	Ketone Resistant	Optional																																															

Nozzle Nuts		
	152290	Standard, zinc plated steel
	152496	Standard, stainless steel
	123200	Anti-foul, stainless steel
	152926	Extended guns

