



1625720
Adapter Injector Kit

Threaded Adapter

Machined Injector

1

-014, Conductive O-Ring

2

Conductive O-Ring



Pump Throat Kit - 10 Pack

1625721
Throat Tivar® Assembly (10 Pack)

1625722
Throat PTFE Assembly (10 Pack)

NOTE: Pumps are shipped with standard Tivar venturi throats, PTFE throats must be ordered separately.

1625719
Throat Holder Kit

1

-014, Conductive O-Ring

Pump Nut

O-Ring Kit 1623816	
1	x 20
2	x 20

O-Ring Reference Chart
(Use this chart to help in the identification of O-rings.)
(Print requirements for proper sizing: 11x17, landscape)

-014, Conductive O-rings
1

Conductive O-rings
2

Optional Tools

1097913
Injector Extraction Tool

152999
Wrench